

IQ NAILS

DRILL BITS FOR NAIL AND FOOT CARE

MADE IN POLAND



Duo-Bits

Golden Line

DOUBLE-GRIT DIAMOND DRILL BIT WITH AN INNOVATIVE NANO LAYER ZRN
WHICH PROTECTS THE DRILL BIT AND EXTENDS ITS LIFETIME.



001.FM.025_ZRN

BALL SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE STANDARD DIAMOND GRIT WITH ADDITIONAL FINE DIAMOND GRIT.



243R.FM.021_ZRN

LONG FLAME WITH ROUND END SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE STANDARD DIAMOND GRIT WITH ADDITIONAL FINE DIAMOND GRIT.



257.FM.045_ZRN

LARGE FLAME SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE STANDARD DIAMOND GRIT WITH ADDITIONAL FINE DIAMOND GRIT.



171.FM.031_ZRN

SHORT FLAT CONE SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE STANDARD DIAMOND GRIT WITH ADDITIONAL FINE DIAMOND GRIT.



172.FM.031_ZRN

LONG FLAT CONE SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE STANDARD DIAMOND GRIT WITH ADDITIONAL FINE DIAMOND GRIT.



001.EFF.023_ZRN

BALL SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE FINE DIAMOND GRIT WITH ADDITIONAL EXTRA FINE DIAMOND GRIT.



243R.EFF.023_ZRN

LONG FLAME WITH ROUND END SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE FINE DIAMOND GRIT WITH ADDITIONAL EXTRA FINE DIAMOND GRIT.



254.EFF.023_ZRN

FLAME SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE FINE DIAMOND GRIT WITH ADDITIONAL EXTRA FINE DIAMOND GRIT.



172.EFF.031_ZRN

LONG FLAT CONE SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE FINE DIAMOND GRIT WITH ADDITIONAL EXTRA FINE DIAMOND GRIT.



260.EFF.060_ZRN

LARGE FLAME SHAPED DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE FINE DIAMOND GRIT WITH ADDITIONAL EXTRA FINE DIAMOND GRIT.

Cool Bits

Premium Line

THE COOL BITS ARE AN EXTREMELY INNOVATIVE AND USEFUL SOLUTION IN E-FILE MANICURE. NON-STANDARD CROSS CUTS ON THE SURFACE OF DIAMOND CUTTERS WITH A DEPTH OF 0.1 MM PROVIDE BETTER AIR CIRCULATION, WHICH PREVENTS OVERHEATING OF THE TREATED SURFACE AND ENSURES BETTER COMFORT DURING THE PROCESS, ESPECIALLY DURING FAST SALON WORK.

ADDITIONALLY, THE CUTTERS ARE COVERED WITH NANO ZRN COATING, WHICH PROTECTS THE DIAMOND FROM CHIPPING AND EXTENDS ITS LIFE DURING DAILY STERILIZATION AND DISINFECTION PROCESSES. PREMIUM LINE COOL BITS INCLUDES THE MOST POPULAR SHAPES IN TWO DIFFERENT GRADES; BLUE GRADE 524 FOR FASTER WORK AND RED GRADE 514 FOR GENTLE AND CONTROLLED WORK.



001.X514.031_ZRN

DIAMOND DRILL BIT WITH THE SHAPE OF BALL ALLOWS YOU TO QUICKLY, ACCURATELY AND NON-INVASIVELY REMOVE RAISED CUTICLES WITHOUT THE NEED TO USE NIPPERS. THE RED RING STANDS FOR DELICATE DIAMOND GRAIN GRADATION.



001.X524.031_ZRN

DIAMOND DRILL BIT WITH THE SHAPE OF BALL ALLOWS YOU TO QUICKLY, ACCURATELY AND NON-INVASIVELY REMOVE RAISED CUTICLES WITHOUT THE NEED TO USE NIPPERS. THE BLUE STRIPE STANDS FOR UNIVERSAL DIAMOND GRAIN GRADATION.



111.X514.031_ZRN

DIAMOND DRILL BIT WITH THE SHAPE OF CYLINDER FOR QUICK AND PRECISE PREPARATION AND MATTIFYING OF THE NAIL PLATE, SMOOTHING THE CUTICLE AND THE SURFACE OF THE NATURAL NAIL PLATE, CREATING A PERFECT C CURVE. RED GRADE 514 FOR GENTLE AND CONTROLLED WORK.



111.X524.031_ZRN

DIAMOND DRILL BIT WITH THE SHAPE OF CYLINDER WILL SAFELY AND GENTLY MOVE AWAY THE CUTICLES AND GROWN EPITHELIA. IT GUARANTEES FAST AND PRECISE WORK. BLUE GRADE 524 FOR FASTER WORK



173.X514.031_ZRN

DIAMOND DRILL BIT WITH THE SHAPE OF A FLAT END CONE FOR PUSHING, MATTIFYING, REMOVAL OF EPITHELIUM, FOR CORRECTING THE SHAPE. RED GRADE 514 FOR GENTLE AND CONTROLLED WORK.



173.X524.031_ZRN

DIAMOND DRILL BIT WITH THE SHAPE OF A FLAT END CONE FOR PUSHING, MATTIFYING, REMOVAL OF EPITHELIUM, FOR CORRECTING THE SHAPE. BLUE GRADE 524 FOR FASTER WORK



243R.X514.023_ZRN

DIAMOND DRILL BIT IN THE SHAPE OF A LONG FLAME FOR: LIFTING THE CUTICLES, REMOVING THE EPITHELIA AND ALIGNING THE EDGES OF THE NAIL. RED GRADE 514 FOR GENTLE AND CONTROLLED WORK.



243R.X524.023_ZRN

DIAMOND DRILL BIT IN THE SHAPE OF A LONG FLAME FOR: LIFTING THE CUTICLES, REMOVING THE EPITHELIA AND ALIGNING THE EDGES OF THE NAIL. BLUE GRADE 524 FOR FASTER WORK



Duo bits




New DUO - BITS SET

THE SET IS A COLLECTION OF A NEW VERSION OF THE MOST POPULAR SHAPES, CREATED FOR INCREASED COMFORT. FOR MAXIMIZED EFFICIENCY AND LIFETIME OF THE BITS, IQNAILS HAS CREATED A DOUBLE-GRIT DIAMOND DRILL BIT, BY FILLING THE SPACES BETWEEN THE STANDARD DIAMOND GRIT WITH ADDITIONAL FINE DIAMOND GRIT. THE SET INCLUDES AN INNOVATIVE CARBIDE DRILL BIT TO WORK IN BOTH DIRECTIONS. ALL DRILL BITS IN THE SET ARE LAYERED WITH AN INNOVATIVE NANO LAYER ZRN WHICH PROTECTS THE DRILL BIT AND EXTENDS ITS LIFETIME.

THE SET CONTAINS:

- FIRE** - CONE SHAPED CARBIDE DRILL BIT WITH CUTS DEDICATED FOR BOTH RIGHT AND LEFT-HANDED STYLISTS,
- LONGER STRONGER** - LONG FLAT CONE SHAPED DIAMOND DRILL BIT,
- HI SHORTY** - SHORT FLAT CONE SHAPED DIAMOND DRILL BIT,
- FAT BOY** - VERY LARGE FLAME SHAPED DIAMOND DRILL BIT,
- LADY FLAME** - FLAME SHAPED DIAMOND DRILL BIT,
- CUTICLE QUEEN** - BALL SHAPED DIAMOND DRILL BIT.




Cool bits



New COOL - BITS SET

THE COOL BITS SET IS AN EXTREMELY INNOVATIVE AND USEFUL SOLUTION IN E-FILE MANICURE. NON-STANDARD CROSS CUTS ON THE SURFACE OF DIAMOND CUTTERS WITH A DEPTH OF 0.1 MM PROVIDE BETTER AIR CIRCULATION, WHICH PREVENTS OVERHEATING OF THE TREATED SURFACE AND ENSURES BETTER COMFORT DURING THE PROCESS, ESPECIALLY DURING FAST SALON WORK.

ADDITIONALLY, THE CUTTERS IN THE SET ARE COVERED WITH NANO ZRN COATING, WHICH PROTECTS THE DIAMOND FROM CHIPPING AND EXTENDS ITS LIFE DURING DAILY STERILIZATION AND DISINFECTION PROCESSES. THE SET INCLUDES THE MOST POPULAR SHAPES IN TWO DIFFERENT GRADES; BLUE GRADE 524 FOR FASTER WORK AND RED GRADE 514 FOR GENTLE AND CONTROLLED WORK.



Black Diamond

Manicure Set
by Marcela

IQ
NAILS



QUICK & EASY SET BY MARCELA DEDICATED TO ENHANCE THE PACE OF WORK

THE QUICK & EASY MANICURE SET BY MARCELA HAS BEEN CREATED IN COLLABORATION WITH MARCELA, A WELL KNOWN POLISH NAILART INSTRUCTOR, TO ENHANCE THE PACE OF WORK. THE SET CONTAINS ALL DRILL BITS NEEDED TO PERFORM: EXTENSION, INFILL, COMBINED MANICURE AND CORRECTION OF DOWNWARD GROWING NAILS.

ALL DRILL BITS IN THE SET HAVE BEEN LAYERED WITH A SPECIAL GOLD ZRN LAYER, WHICH GUARANTEES:

- MINIMIZED OVERHEATING EFFECT,
- BETTER BIOCOMPATIBILITY – RECOMMENDED FOR PEOPLE ALLERGIC TO METALS,
- EXTENDED DRILL BIT LIFETIME, MORE ECONOMICAL OPERATION, GREATER RESISTANCE TO STERILIZATION SIDE EFFECTS AND QUICK REMOVAL OF DECORATIONS SUCH AS GLITTER AND CRYSTALS.



QUICK & EASY MANICURE BY MARCELA

THE SET CONTAINS TWO CARBIDE DRILL BITS FOR MATERIAL REMOVAL AND FOUR DIAMOND DRILL BITS FOR CUTICLES.

IQN SET MANI 27



Pink Diamond



ALL IN ONE SET BY MATKA POLKA HYBRYDOWA

IN THE SALON, IT IS NOT ONLY THE PACE OF WORK THAT COUNTS, BUT MOST OF ALL THE EFFECT. PERFECT CUTICLES ARE IMPOSSIBLE TO ACHIEVE WITHOUT THE HIGHEST QUALITY TOOLS. A WELL-THOUGHT-OUT AND TESTED FOR YEARS SET OF CUTTERS BY MATKA POLKA GEL POLISH ALLOWS YOU TO WORK WITH ANY TYPE OF CUTICLE AND WITH ANY MATERIAL DURING MANICURE AND PEDICURE.

ALL DIAMOND AND CARBIDE DRILL BITS IN THE SET HAVE BEEN LAYERED WITH A SPECIAL GOLD ZRN LAYER, WHICH GUARANTEES:

- MINIMIZED OVERHEATING EFFECT,
- BETTER BIOCOMPATIBILITY – RECOMMENDED FOR PEOPLE ALLERGIC TO METALS,
- EXTENDED DRILL BIT LIFETIME, MORE ECONOMICAL OPERATION, GREATER RESISTANCE TO STERILIZATION SIDE EFFECTS AND QUICK REMOVAL OF DECORATIONS SUCH AS GLITTER AND CRYSTALS.



ALL IN ONE SET BY MATKA POLKA HYBRYDOWA

THE SET CONTAINS ONE CARBIDE DRILL BIT FOR MATERIAL REMOVAL, FOUR DIAMOND AND ONE STONE DRILL BITS FOR CUTICLES.

IQN SET MANI 32



EXPERT FAMILY



EXPERT DRILL BIT SETS LINE

EXPERT SETS ARE A COMPREHENSIVE COMBINATION OF DRILL BITS FOR CUTICLE CARE AND QUICK REMOVAL OF GEL POLISH, GEL AND ACRYLIC, CREATED FOR DEMANDING STYLISTS.

THE SETS ARE DESIGNED FOR BOTH BEGINNER AND ADVANCED NAIL STYLISTS. THE EXPERT SET LINE INCLUDES DRILL BITS WITH SPECIAL NANO LAYERS, WHICH GUARANTEE:

- MINIMIZED OVERHEATING EFFECT (CRUCIAL FEATURE WHEN PREPARING CUTICLES AND WORKING ON THE NATURAL NAIL PLATE, FASTER MATERIAL REMOVAL),
- REDUCED RISK OF MATERIAL STICKING TO THE DRILL BITS,
- THREE TIMES LONGER LIFETIME OF THE DRILL BIT,
- COST-EFFECTIVE SOLUTION,
- GREATER RESISTANCE TO STERILIZATION SIDE EFFECTS.
- QUICK REMOVAL OF NAIL DECORATIONS SUCH AS GLITTER OR CRYSTALS.



EXPERT I

THE SET CONTAINS A CARBIDE DRILL BIT FOR MATERIAL REMOVAL, FOUR DIAMOND DRILL BITS FOR CUTICLES, INCLUDING TWO WITH A CHAMELEON NANO LAYER AND ANOTHER TWO WITH A ZIRCONIUM NANO LAYER.

IQN SET MANI 11



266.190.060



173.524.016_ZRN



243R.524.021_KVN



257.514.018_KVN



001.524.021_ZRN

EXPERT II

THE SET CONTAINS A CARBIDE DRILL BIT WITH A ZIRCONIUM NANO LAYER FOR REMOVING MATERIALS AND FOUR DIAMOND CUTICLE DRILL BITS, INCLUDING TWO WITH A CHAMELEON NANO LAYER. THIS SET IS PERFECT FOR STYLISTS FAMILIAR WITH THE DRY MANICURE TECHNIQUE AND WHO WANT TO START A NEW EXPERIENCE WITH MORE ADVANCED MANICURE TECHNIQUES.

IQN SET MANI 12



EXPERT GOLD POWER

THE SET CONTAINS A CARBIDE DRILL BIT FOR MATERIAL REMOVAL AND FIVE DIAMOND CUTICLE DRILL BITS. ALL OF THEM ARE LAYERED WITH A ZIRCONIUM NANO LAYER. THIS SET HAS BEEN CREATED FOR STYLISTS WHO APPRECIATE FAST WORK.

IQN SET MANI 14



EXPERT PRO

THE SET CONTAINS A CARBIDE DRILL BIT FOR MATERIAL REMOVAL AND FIVE DIAMOND CUTICLE DRILL BITS, INCLUDING THREE WITH A CHAMELEON NANO LAYER AND TWO WITH A NANO ZIRCONIUM LAYER. THIS SET HAS BEEN COMPOSED TO EXPEDITE WORKING TIME EFFECTIVELY.

IQN SET MANI 13



EXPERT BASIC

THE SET CONTAINS A CARBIDE DRILL BIT FOR MATERIAL REMOVAL AND FOUR DIAMOND CUTICLE DRILL BITS, WHICH HAVE BEEN LAYERED WITH AN INNOVATIVE SAPPHIRE NANO LAYER. IN ADDITION, THE CUTICLES WILL BE TAKEN CARE OF BY THE CERAMIC CUTTER INCLUDED IN THE SET. THIS SET IS ESPECIALLY RECOMMENDED FOR BEGINNERS, LOVERS OF CONTROLLED AND SAFE WORK.

IQN SET MANI 15



Paulina Pastuszek


DRILL BIT SETS BY PAULINA PASTUSZAK AN EXPERT IN THE FIELD OF NAIL STYLING



PAULINA PASTUSZAK, THE WINNER OF THE "NAIL ART WORLD CUP" IN MUNICH AND A SEVEN-TIME FINALIST OF THE POLISH CHAMPIONSHIPS, ALMOST 30-TIME THE JURY IN POLISH AND INTERNATIONAL NAIL STYLING CHAMPIONSHIPS, AUTHOR OF OVER 400 ARTICLES IN THE TRADE PRESS AROUND THE WORLD,

THE COLLABORATION WITH PAULINA RESULTED IN THE CREATION OF SPECIAL DRILL BIT SETS, WHICH WILL FACILITATE THE EVERYDAY WORK OF STYLISTS. THEY PRESENT A COMPLETELY NEW LOOK ON THE DRY MANICURE!

THREE UNIQUE SETS WERE CREATED TOGETHER WITH PAULINA, WITH GREAT ATTENTIONS TO EVERY DETAIL WHEN COMPOSING THE CONTENT. AMONG THE SETS, YOU WILL FIND AN INNOVATIVE DRILL BIT, CREATED BY OUR TEAM OF EXPERTS ON THE BASIS OF PAULINA'S ORIGINAL DESIGN.



KMN-147XUF HP.060_PP

THE DRILL BIT WAS SUBJECTED TO THE PATENT PROCEDURE DUE TO AN UNUSUAL SOLUTION PROTECTING LATERAL AND PROXIMAL FOLDS. THE CYLINDER-SHAPED CARBIDE DRILL BIT WITH A SAFE TOP AND AN INNOVATIVE SAFE END AT THE BOTTOM OF THE WORKING PART ENSURES SAFE AND NON-INVASIVE OPERATION. THE DRILL BIT HAS A SPECIAL NANO KMN COATING.

SENSITIVE

SOME OF YOUR CLIENTS' CUTICLES ARE VERY THIN, TOO FLEXIBLE, PULLED UP OR TORN? THIS CASE REQUIRES THE USE OF DRILL BITS THAT WILL ALLOW YOU TO WORK DECREASING THE RISK OF BREAKING THE CONTINUITY OF THE TISSUE. THE DRILL BITS INCLUDED IN THIS SET WILL TAKE CARE OF THE PERFECT MANICURE FOR DELICATE HANDS WITH GREAT ATTENTION.

IQN SET MANI 20



KMN-147XUF-HP.060_PP



KMN-274XF-HP.014



111.514.025_KMN



190NC.071.012_KMN



001.524.060_KMN

MAGIC

CARBIDE DRILL BITS HAVE NO MERCY ON EPITHELIUM, CALLUSES, AND IMPERFECTIONS AT THE PROXIMAL FOLD, WHILE GENERATING NOTICEABLY LESS DUST THAN DIAMOND DRILL BITS. THANKS TO THIS, YOUR WORK WILL BE QUICK AND EFFICIENT.

IQN SET MANI 21




KMN-275QSG-HP.060



KMN-225XEF-HP.023



KMN-198XF-HP.023



001.524.029_KMN



260.504.060_KMN


DEEP GEL MANICURE

THIS SET WILL PERFECTLY SATISFY YOUR NEED FOR A SLIGHTLY DEEPER INTERVENTION INTO THE PERIUNGUAL SHAFT AREA USING A DRILL MACHINE. YOU WILL GET RID OF IMPERFECTIONS QUICKLY, EFFICIENTLY AND NON-INVASIVELY - WITHOUT HAVING TO USE CUTTING TOOLS. WE RECOMMEND THIS SET FOR ADVANCED STYLISTS.


IQN SET MANI 19



KMN-147XUF-HP.060_PP



KMN-225X-HP.016



KMN-001-HP.021



001.524.012_KMN



KMN-198XS-HP.016



260.524.060_KMN




111.X514.050_ZRN



111.X524.050_ZRN



111.X514.031_ZRN



111.X524.031_ZRN

DRILL BIT LIKE MANDREL STANDARD AND MINI

THIS CYLINDER - SHAPED BIT WITH SPECIAL CUTS AND ADDITIONAL FINE DIAMOND GRIT IS PERFECT FOR PERFORMING SURFACE WORK, PREPPING AND SMOOTHING THE NAIL BED ON EACH OF YOUR CLIENTS! IF YOU PREFER ECO SOLUTION THIS BIT IS MADE FOR YOU.


Final Touch Set

PERFECT CUTICLE ROUTINE



FINAL TOUCH CUTICLE OIL 30 ML

REMEMBER, STYLISTS WON'T BE ABLE TO ACHIEVE A PERFECT NAIL LOOK WITHOUT HAVING PERFECTLY SMOOTHED CUTICLES. FINAL TOUCH CUTICLE OIL IS A COMPREHENSIVE COMBINATION THAT CONTAINS, AMONG OTHERS, VALUABLE OILS: GRAPE SEED AND PEACH SEED OILS, AS WELL AS CAREFULLY SELECTED EMOLLIENTS THAT SMOOTHENS THE SKIN AFTER A MANICURE SESSION. IT IS BEST TO USE OLIVE OIL DURING THE FINAL POLISHING OF THE CUTICLES USING THE "WET" METHOD FOR ACHIEVING PERFECTLY SMOOTHED CUTICLES. THE SELECTION OF FRAGRANCE IN THE OIL, GIVES THE VIBES OF EXOTIC BEACH, SAND, SEASIDE SUN AND SWEET FRUIT, WHICH WILL ADDITIONALLY EMBRACE YOUR HANDS WITH A BEAUTIFUL EXOTIC SCENT TILL THE END OF THE DAY. THE USE OF OLIVE OIL AT THIS STAGE OF NAIL PREPARATION FOR THE FINAL STYLING ELIMINATES THE NEED TO LUBRICATE THE CUTICLES AT THE VERY END OF THE MANICURE SESSION.



HOW TO USE: MASSAGE THE OIL ONTO PRE-PREPARED CUTICLES AND NAIL FOLDS. USE THE OIL SIMULTANEOUSLY WITH A FINE DRILL BIT FOR POLISHING. ALLOW IT TO BE ABSORBED FOR A BETTER MOISTURIZING EFFECT. IF YOU ARE GOING TO PERFORM GEL POLISH STYLING, DEGREASE THE NAIL THOROUGHLY.

SOMETHING FOR *left-handed*



LEFT-HANDED STYLISTS SETS

THE SETS INCLUDE DRILL BITS WITH SPECIAL CUTS DESIGNED FOR LEFT-HANDED STYLISTS. THE SETS INCLUDE CARBIDE AND CERAMIC DRILL BITS - WITH SPECIAL CUTS FOR WORKING LEFTWARD AND DIAMOND DRILL BITS FOR CUTICLES.

MINI LEFT-HANDED STYLIST SET

THE SET INCLUDES: ONE CARBIDE AND ONE CERAMIC DRILL BITS WITH SPECIAL CUTS FOR WORKING LEFTWARD AND DIAMOND DRILL BITS FOR CUTICLE PREPERATION.

IQN SET MANI 23



LEFT-HANDED STYLIST SET

THE SET INCLUDES TWO CARBIDE DRILL BITS WITH SPECIAL LEFT CUTTING FLUTES AND DIAMOND DRILL BITS FOR CUTICLE PREPARATION.


IQN SET MANI 24



Gold SET



PREMIUM "GOLD" DRILL BIT SETS FOR PROFESSIONALS



GOLD MAXI

UNIVERSAL DRILL BIT SET FOR STYLING REMOVAL - GEL POLISH/ GEL/CUTICLES. THE SET INCLUDES 2 CARBIDE DRILL BITS AND 4 DIAMOND DRILL BITS FOR PROCESSING BOTH DELICATE AND MORE CALLOUSED CUTICLES.

IQN SET MANI 17



ZRN-201XF HP.060



ZRN-274X HP.050




001.514.018_ZRN



243R.514.021_ZRN



001.524.018_ZRN




243R.524.021_ZRN

GOLD SET GEL POLISH/ GEL

UNIVERSAL DRILL BIT SET FOR GEL POLISH AND GEL REMOVAL. THE SET INCLUDES 4 CARBIDE DRILL BITS THAT ALLOWS REMOVAL OF BOTH SOFT AND MEDIUM-SOFT MATERIALS.

IQN SET MANI 18



ZRN-201XF HP.060




ZRN-201X HP.040



ZRN-274X HP.050



ZRN-274V HP.050



GOLD MINI

UNIVERSAL DRILL BIT SET FOR STYLING REMOVAL - GEL POLISH/GEL/CUTICLES. THE SET INCLUDES 2 CARBIDE DRILL BITS AND 2 DIAMOND DRILL BITS FOR CUTICLE PREPARATION.

IQN SET MANI 16




ZRN-201XF HP.060



ZRN-274X HP.050



001.524.018_ZRN




243R.524.021_ZRN

UNIVERSAL



UNIVERSAL SETS

CHOOSING DRILL BIT SETS IS A COST-EFFECTIVE CHOICE THAT INCLUDES ALL THE NECESSARY ITEMS IN ONE PACK. WITH THESE SETS YOU WILL COMPREHENSIVELY PREPARE CUTICLES AND REMOVE ALL MATERIALS - GEL POLISH, GEL, ACRYLIC AND DIP POWDER.



Manicure Set

7 GEL POLISH DRILL BITS / GEL / CUTICLE

THE GEL POLISH, GEL AND CUTICLE SET INCLUDES A CARBIDE DRILL BIT WITH A SAFE WORKING UPPER PART AND 6 DIAMOND DRILL BITS FOR CUTICLE REMOVING.


IQN SET MANI 5



UNIVERSAL STYLING REMOVAL SET

UNIVERSAL SET OF DRILL BITS FOR REMOVING STYLING - GEL POLISH, GEL, ACRYLIC AND DIP POWDER. THE SET INCLUDES: 3 CARBIDE AND 1 CERAMIC DRILL BITS.

IQN SET MANI 8





PEDI cure



bezpieczny zestaw
rezów dla
diabetyków

IQ
NAILS

rekomenduje


Milena Skalska



MEDICAL DRILL BITS FOR PODIATRISTS

IQNAILS COLLABORATES NOT ONLY WITH EXPERTS IN THE FIELD OF NAIL STYLING, BUT IN PODOLOGY AS WELL. IN OUR PODOLOGICAL OFFER YOU WILL FIND DRILL BIT SETS CREATED IN COLLABORATION WITH PODOLAND, MILENA SKALSKA AND DANIEL ABRATKIEWICZ.


PROFESSIONAL DRILL BITS FOR PODIATRISTS WILL ENSURE THE SAFETY OF CLIENTS AND FACILITATE YOUR WORK ON NAILS THAT HAVE UNDERGONE VARIOUS CHANGES.



"MUST-HAVE-SET" FOR PODIATRISTS

A COMPREHENSIVE SET OF 13 MEDICAL DRILL BITS RECOMMENDED BY PODOLAND. APPLICATION: WORK RECONSTRUCTING ACRYLIC NAILS, REMOVAL OF CORNEA, CALLUSES AND CORNS, CLEANING OF GREEN AND BLUE BACTERIA, CLEANING NAIL FOLD, CLEANING SUBUNGUAL SPACES, TREATMENT SUBUNGUAL HEMATOMAS, TREATMENT OF CRACKS IN THE HEELS, WORKING ON AND MATTIFYING THE NAIL PLATE, NAIL SHORTENING, PARENTHESSES REMOVAL ERADITION ONYCHOLYSIS AND MYCOSIS.

IQN SET PODO 2





„BECAUSE NOT EVERYTHING CAN BE PAINTED OVER ” BY DANIEL ABRATKIEWICZ
IF YOU'VE EVER WONDERED WHICH SPECIFIC DRILL BIT TO USE WHEN A FEMALE OR MALE CLIENT WITH A DISEASED NAIL CAME TO YOU – THIS TIME YOU WON'T HAVE TO! ASSUMING THAT WE KNOW EXACTLY HOW TO HELP THE CLIENT AND WE DECIDE TO DEAL WITH THEIR PROBLEM, THIS SET IS DESIGNED FOR YOU, STYLISTS AND PODIATRISTS. IT WILL HELP YOU DEAL WITH ALL ADVERSITIES. HERE YOU WILL GET EVERYTHING NEEDED.

IQN SET PODO 2



157.DA.524.050



199.524.012.DLC



166.524.018.DLC



DLC-001 HP.012



001.524.018.DLC



HEALTHY FEET & HANDS

COMPREHENSIVE SET OF 6 MEDICAL DRILL BITS RECOMMENDED BY PODOLAND. FOR CLEANSING: MECHANICALLY DAMAGED NAILS, DISEASED NAILS, PSEUDOMONAS, GREEN BACTERIA, INGROWN NAILS, CRACKING HEELS, REMOVAL OF CORNS AND CALLUSES, AND FOR REMOVING STYLINGS AND DEVELOPING CUTICLES.

IQN SET PODO 3



227.4XF HP 060_P



110.534.050



199.524.012




001.524.018



257.534.023



001.175.012



SET OF SAFE DRILL BITS FOR DIABETICS BY MILENA SKALSKA

PROBLEMS WITH WOUNDS HEALING AND IMPAIRED SENSATION ARE COMMON IN PATIENTS WITH DIABETES. HAVING SUCH PATIENTS IN MIND, A SPECIALIZED AND SAFE SET OF DRILL BITS WAS CREATED IN COOPERATION WITH AN EXPERT IN THE FIELD OF PODOLOGY – MILENA SKALSKA.

IQN SET PODO 4



001.514.021_ZRN



001.514.031_ZRN



157TR.514.050_ZRN



157SR.514.050_ZRN



140.514.023_ZRN



001.514.060_ZRN

Innovative bits by Daniel Abratkiewicz



4 in 1



AN INNOVATIVE UNIQUE DRILL BIT WITH A DIAMOND TIP, THAT ALLOWS FAST AND SAFE NAIL SHORTENING WITHOUT THE NEED TO REPLACE THE DRILL BIT FOR THIS PURPOSE.

FINE CROSSCUTS SIMILAR TO THE CERAMIC DRILL BIT FROM „BECAUSE NOT EVERYTHING CAN BE PAINTED OVER" SET, ALLOWS FAST PROCESSING OF HARD MATERIALS AS ACRYLIC, GEL AND DENTAL COMPOSITE.

ITS CYLINDRICAL SHAPE ALLOWS TO ACCESS HARD-TO-REACH AREAS OF THE NAIL PLATE. IN ADDITION TO THE EXCEPTIONAL DURABILITY THAT ALLOWS TO REMOVE MATERIALS THAT ARE TOO HARD FOR OTHER CUTTERS.



THE PODIATRIC DRILL BIT WAS CREATED IN COLLABORATION WITH ZNANY PODOLOG.

THIS BIT HAS FINE SPIRAL CROSSCUTS AS WELL AS A SPECIAL ANGE AND IS USED TO REMOVE DENTAL AND PODIATRY COMPOSITE FROM THE TOENAILS. AND EFFICIENT WORK IN PODIATRY COMPOSITE.

THE TOOL IS MANUFACTURED ENTIRELY FROM FULL CARBIDE WHICH PROVIDES SAFE AND CONTROLLED WORK.

ADDITIONALLY, THE SHARPNESS AND DEPTH OF THE CROSSCUTS ENABLE FAST AND PRECISE WORK WITH THE COMPOSITE AS WELL AS LESIONS ON THE NAIL PLATE.

THE BIT HAS BEEN PROTECTED BY SPECIAL DLC „DIAMOND LIKE CARBON" COATING WHICH IS EXTREMELY HARD AND RESISTANT TO CORROSION AND USAGE.



A DRILL BIT FOR PROCESSING THICKENED LESIONS ON THE FOOT SOLES.

THE DRILL BIT WITH DIAMOND POWDER ENSURES PRECISE REMOVAL OF THICKENINGS, CORNS AND CALLUS. FINE POWDER SMOOTHENS SKIN AFTER THE INTERVENTION OF A SCALPEL OR MORE AGGRESSIVE ROTATING TOOLS DRILL BITS.



V104.144S.514.070



V104.144S.524.070



V104.144S.534.070



V104.144S.544.070



DRILL BIT FOR PROCESSING FOOT SOLES.

THE DRILL BIT WITH DIAMOND POWDER ENSURES PROCESSING THE SURFACE OF FOOT SOLES, REMOVING THICKENINGS, CORNS AND CALLUS. FINE DIAMOND POWDER SMOOTHENS SKIN AFTER THE INTERVENTION OF A SCALPEL OR MORE AGGRESSIVE ROTATING TOOLS DRILL BITS.



V104.144S.514.100



V104.144S.524.100



V104.144S.534.100



V104.144S.544.100



DRILL BIT FOR PROCESSING FOOT SOLES.

THE DRILL BIT WITH DIAMOND POWDER ENSURES PROCESSING THE SURFACE OF FOOT SOLES, REMOVING THICKENINGS, CORNS AND CALLUS. FINE DIAMOND POWDER SMOOTHENS SKIN AFTER THE INTERVENTION OF A SCALPEL OR MORE AGGRESSIVE ROTATING TOOLS DRILL BITS.



V104.144S.514.130



V104.144S.524.130



V104.144S.534.130



V104.144S.544.130

Solo Bestseller



Soft Heels



PodopACKS

drill bits



PACK OF 5 DIAMOND PDIATRY DRILL BITS


OUR OFFER INCLUDES PDIATRY DRILL BIT PACKS WITH THE MOST POPULAR DRILL BITS USED IN PDIATRY LIKE: CYLINDER, NEEDLE, ROUND. THE SETS CONTAIN TWO DIFFERENT VERSIONS : BASIC AND PREMIUM



PODIATRY PACK I

PODIATRY DRILL BIT PACK INCLUDES 3 OF THE MOST POPULAR DRILL BITS USED IN PDIATRY. IN THE SET YOU WILL FIND: CARBIDE ROUND DRILL BIT, DIAMOND ROUND DRILL BIT AND TREPHINE DRILL BIT.

IQN PACK 3



PODIATRY PACK II

PACK OF PODOLOGY DRILL BITS II INCLUDES 3 OF THE MOST POPULAR DRILL BITS USED IN PDIATRY. IN THE SET YOU WILL FIND: CERAMIC PINK DRILL BIT, LARGE CYLINDRICAL DIAMOND DRILL BIT AND DIAMOND DRILL BIT ROUNDED CONE.

IQN PACK 4



PACKS

drill bits



HOME SPA GOLD

DRILL BITS PACK HOME SPA GOLD CONTAINS THREE DRILL BITS WITH GOLDEN ZRN COATING, WHICH GUARANTEES: MINIMIZING THE OVERHEATING EFFECT OF THE NAIL SURFACE, BETTER BIOCOMPATIBILITY OF THE DRILL BIT – RECOMMENDED FOR PEOPLE WHO ARE ALLERGIC TO METALS, LONGER SERVICE LIFETIME OF THE DRILL BIT, MORE ECONOMICAL WORK, HIGHER RESISTANCE TO ALL STERILANTS SIDE EFFECTS.

IQN PACK 2



198.524.023_ZRN



ZRN-147XUF HP.060



001.524.023_ZRN



BASIC HOME

BASIC HOME PACK IS A MUST HAVE TO START WORKING WITH A DRILL MACHINE. SPECIALLY SELECTED DRILL BITS WILL ALLOW YOU TO REMOVE THE OLD STYLING (GEL POLISH/GEL) AND GENTLY REMOVE THE CUTICLES.

IQN PACK 1



201.134.060



198.514.023




001.514.023



RAINBOW

PACK OF THREE COLORED DRILL BITS WITH COATING THAT ENABLES YOU TO DO FULL MANICURE PROCESS.

IQN PACK 18



171.FM.031_ZRN



KMN-137XUF HP.060



243R.514.021

PACKS

drill bits



DRILL BIT PACK FOR LEFT-HANDED STYLISTS
THE PACK FOR LEFT-HANDED STYLISTS INCLUDES: A CARBIDE DRILL BIT WITH SPECIAL CUTS FOR WORKING LEFTWARD AND TWO DIAMOND DRILL BITS FOR PREPARING CUTICLES. THIS PACKAGE IS IDEAL FOR NAIL STYLING AND PODIATRY.

IQN PACK 5



257.524.016



L-274X HP 060



172.524.025



SOFT

WE PRESENT YOU THE SOFT DRILL BITS PACK, CREATED FOR BEGINNER STYLISTS AND PEOPLE WHO LIKE SAFE AND CONTROLLED WORK.

IQN PACK 19



198.514.021



274.110.060



243R.514.023



PACKS OF 5 DIAMOND DRILL BITS

OUR OFFER INCLUDES PACKS OF DIAMOND DRILL BITS PACKED IN 5 PIECES: PEAR, OLIVE, SMALL FLAME, LONG FLAME, TRUNCATED CONE, ROUNDED CONE AND BALL. THE PACKS ARE PERFECT FOR STYLISTS WHO VALUE FAST AND SAFE WORK IN THE SALON. THE PACKAGES HAVE 3 DIFFERENT PACKAGING DESIGNS.



Tool off

WE INTRODUCE BITS THAT WILL BE ABLE TO REPLACE EVERY REGULAR TOOL FOR MANICURE. SAVE YOUR TIME AND MONEY USING DIAMOND BITS FOR SPECIAL TASKS. YOU MIGHT NEED TO PAY A FEW TIMES MORE FOR HIGH QUALITY TOOLS FOR MANICURE AS COMPARED TO ONE BIT. THAT'S WHY WE ARE BRINGING TO YOU REPLACEMENTS FOR PUSHERS, BUFFERS, NIPPERS OR EVEN FILES.



FLAT END TAPER BIT - PUSHER TYPE 173.524.025
THE SIZE OF THE BIT HAS BEEN INCREASED AS COMPARED TO REGULAR TAPER TO HELP YOU WORK FASTER. IT EFFICIENTLY PUSHES UP AND CLEANS THE PLATE EVEN FROM THE ROUGHEST CUTICLES. THE BIT IS RECOMMENDED FOR LONG OR WIDE NAILS.



ROUND END TAPER BIT - BUFFER TYPE 260.514.060
THE BIT SIZE WILL GIVE YOU FAST MATTIFYING OF THE ENTIRE NAIL SURFACE AND THE FINE GRIT WILL ONLY REMOVE EXTERNAL LAYER OF THE NAIL PLATE. THIS SHAPE WILL ALSO BE PERFECT FOR POLISHING CUTICLES.



FLAME SHAPED BIT - UNDER THE SKIN TYPE - 540.514.010
FLAME SHAPED DIAMOND DRILL BIT FOR: LIFTING THE CUTICLES, REMOVING EPIPHYSES, IDEALLY SUITED AS A "CLEANER", WILL REMOVE EPIPHYLES EVEN FROM HARD-TO-REACH PLACES AND ALLOW YOU TO PUT GEL POLISH DEEPER.



CYLINDER BIT - MANDREL TYPE - 113.X514/524.060
THIS CYLINDER - SHAPED BIT WITH SPECIAL CUTS AND ADDITIONAL FINE DIAMOND GRIT IS PERFECT FOR PERFORMING SURFACE WORK, PREPPING AND SMOOTHING THE NAIL BED ON EACH OF YOUR CLIENTS! IF YOU PREFER ECO SOLUTION THIS BIT IS MADE FOR YOU.



ROUND BIT - NIPPER/SCISSORS TYPE - 001.524.031
IF YOU DON'T WANT TO USE CUTTING TOOLS OR YOUR CLIENTS ARE CONCERNED ABOUT CUTICLES GROWING BACK FASTER, WE HAVE A PERFECT SOLUTION FOR YOU - OUR ROUND BIT. IT WILL REMOVE EXTRA CUTICLES WITHOUT LEAVING UNEVEN SURFACES OR CUTS.

File off

TAPER AND CYLINDER BITS - FILE TYPE - OUR BITS WILL PROVIDE YOU WITH BOTH HIGH QUALITY AND EFFICIENT SOLUTION REPLACING THE WORK DONE SOMETIMES BY MANY FILES. MULTIPLE DENSE CUTS ON OUR BITS WILL ENSURE SMOOTH SURFACE AND DIFFERENT SHAPES AND ENDS WILL GIVE YOU THE FLEXIBILITY AND COMFORT OF WORK. YOU CAN USE THE BITS FOR GEL POLISH, GEL AND ACRYLIC MANICURE.




sofa bestseller



201.110.060

CARBIDE DRILL BIT WITH EXTRA FINE CROSS CUTS. FOR REMOVAL OF GEL POLISH OFF THE NATURAL NAIL PLATE, AS WELL AS GEL AND ACRYLIC SURFACES. IT ALSO CORRECTS THE GEL AND ACRYLIC AROUND THE NAIL FOLDS.




201.190.040/060

CARBIDE DRILL BIT IN THE SHAPE OF A LONG CONE WITH STANDARD CROSS CUTS FOR SAFE REMOVAL OF GEL POLISH AND GEL.



274.140.050/060

CARBIDE DRILL BIT WITH FINE CROSS CUTS. FOR REMOVAL OF GEL POLISH OFF THE NATURAL NAIL PLATE, AS WELL AS GEL AND ACRYLIC SURFACES. IT ALSO CORRECTS THE GEL AND ACRYLIC AROUND THE NAIL FOLDS.




KMN-275QSG HP.060

DRILL BIT FOR QUICK AND EASY REMOVAL OF DIP POWDER MANICURE. A SPECIAL TYPE OF CUTS MAKES IT VERY SHARP AND THE COLORED KMN COATING MINIMIZES THE EFFECT OF OVERHEATING DURING WORK.



ZRN-266XUF HP.060

THE CHRISTMAS TREE-SHAPED CUTTER WITH CROSS-CUTS, WHICH WILL ENSURE THE SMOOTHNESS OF THE SURFACE, WILL ALSO BE AN IDEAL SOLUTION FOR BEGINNERS. THE UPPER PART OF THE CUTTER HAS VERY SMALL CUTS DESIGNED FOR SAFE WORK BY THE CUTICLES. DUE TO ITS DELICACY, IT WILL BE PERFECT FOR REMOVING GEL MANICURE.



ZRN-147XUF HP.060

A CYLINDER-SHAPED CARBIDE DRILL BIT WITH A NON-WORKING UPPER PART THAT WILL: REMOVE THE OLD STYLING, REPLACE THE FILE DURING GEL SHAPING, SHORTEN THE NATURAL AND EXTENDED NAIL AND ENSURE SAFETY WITHIN THE CUTICLES.



ZRN-221XUF HP.060

THE CONE-SHAPED CARBIDE DRILL BIT WITH A SAFE TOP THAT GUARANTEES SAFE WORK IN THE AREA OF THE NAIL FOLDS, EFFECTIVELY AND QUICKLY COPES WITH THE HYBRID AND GEL ON A NATURAL PLATE, HAS BEEN COVERED WITH A GOLD NANO-COATING.



ZRN-274I HP.060

CARBIDE DRILL BITS IN THE SHAPE OF A CONE, WHICH: WILL REMOVE THE OLD STYLING, LIFT THE NAILS GROWING DOWNWARDS AND SHORTEN THE NATURAL AND EXTENDED NAIL.



KMN-292X HP.040

CARBIDE DRILL BIT WITH STANDARD CROSS CUTS. FOR REMOVAL OF GEL POLISH AND GEL OFF THE NATURAL NAIL PLATE, AS WELL AS GEL AND ACRYLIC SURFACES. IT ALSO CORRECTS THE GEL AND ACRYLIC AROUND THE NAIL FOLDS.



266.190.060

CARBIDE DRILL BIT IN THE SHAPE OF CHRISTMAS TREE WITH STANDARD CROSS CUTS FOR SAFE REMOVAL OF GEL POLISH AND GEL.

sofa bestseller



ZRN-201RLF HP.060

CARBIDE DRILL BIT WITH A ROUNDED CONE - ITS SHAPE WILL ALLOW YOU TO SAFELY PROCESS EVEN THE SMALLEST TOENAILS. SIMPLE, DENSE CUTS WILL EFFECTIVELY AND QUICKLY COPE WITH GEL POLISH, GEL AND ACRYLIC. ZRN GOLD COATING ENSURES THAT THE OVERHEATING EFFECT OF THE PROCESSED SURFACE IS MINIMIZED.



KMN-274V HP.050/060

CHAMELEON EFFECT DRILL BIT WITH CROSSWISE CUTS FOR GENTLE REMOVAL OF GEL POLISH AND SMOOTHING THE GEL.



ZRN-243R XEF HP.018

FLAME SHAPED CARBIDE DRILL BIT WITH A SAFE TOP FOR LIFTING CUTICLES, IT GENTLY REMOVES EPITHELIA. THE GOLD COATING ZRN MINIMIZES THE EFFECT OF OVERHEATING DURING WORK.



113.524.060

PODIATRY DIAMOND BIT FOR CLEANSING AND SMOOTHING PROXIMAL PLAQUE AND CUTICLES. IT LEAVES THE NAIL PLATE SMOOTH, TO BE POLISHED. THE BLUE RING STANDS FOR UNIVERSAL DIAMOND GRAIN GRADATION.



113.514.060

DIAMOND DRILL BIT FOR MATTIFYING NATURAL NAIL PLATE AND CUTICLES. LEAVES THE NAIL PLATE SMOOTH. RED COLORED RING INDICATES FINE GRADATION OF THE DIAMOND GRIT.



201.140.040/060

CARBIDE DRILL BIT WITH FINE CROSS CUTS FOR GEL AND ACRYLIC CORRECTION, GENTLE REMOVAL OF GEL POLISH; GENTLE TREATMENT OF CALLOUSED NAILS, THICK SKIN AND CORNS; TREATMENT OF NAIL FUNGUS.



KMN-274X HP.050/060

CHAMELEON EFFECT DRILL BIT WITH STANDARD CROSS CUTS FOR QUICK REMOVAL OF GEL POLISH, GEL AND ACRYLIC; DEVELOPED CORNS, THICKENINGS ON THE SKIN, NAILS WITH MYCOSIS.



141.514.023

LONG ROUNDED CYLINDER FOR QUICK AND PRECISE PREPARATION AND MATTIFYING OF THE NAIL PLATE, SMOOTHING THE CUTICLE AND THE SURFACE OF THE NATURAL NAIL PLATE, CREATING A PERFECT C CURVE.



260.514.045-060

DIAMOND DRILL BIT WITH THE SHAPE OF A LARGE FLAME WITH FOR REMOVING CUTICLES, WHICH REPLACES SCISSORS OR CLIPPERS.



238.514/524.023/031

DIAMOND DRILL BIT ALLOWS YOU TO QUICKLY, ACCURATELY AND NON-INVASIVELY REMOVE RAISED CUTICLES WITHOUT THE NEED TO USE NIPPERS.

sofa bestseller



001.514.010-060

THIS BIT ALLOWS YOU TO QUICKLY, ACCURATELY AND NON-INVASIVELY REMOVE RAISED CUTICLES WITHOUT THE NEED TO USE NIPPERS. THE RED RING STANDS FOR DELICATE DIAMOND GRAIN GRADATION.



001.524.010-060

THIS BIT ALLOWS YOU TO QUICKLY, ACCURATELY AND NON-INVASIVELY REMOVE RAISED CUTICLES WITHOUT THE NEED TO USE NIPPERS. THE BLUE STRIPE STANDS FOR UNIVERSAL DIAMOND GRAIN GRADATION.



140.524.023

ROUNDED CYLINDER FOR QUICK AND PRECISE PREPARATION AND MATTIFYING OF THE NAIL PLATE, SMOOTHING THE CUTICLE AND THE SURFACE OF THE NATURAL NAIL PLATE, CREATING A PERFECT C CURVE.



198.524.012/014/018/023/040

LONG NEEDLE-SHAPED DIAMOND DRILL BIT FOR REMOVING EPITHELIUM FROM HARD-TO-REACH PLACES. THE BIT IS PERFECT FOR SAFELY LIFTING VERY DELICATE CUTICLES. IT IS INDISPENSABLE IN THE "PAINTING UNDER THE SKINS" TECHNIQUE.



199.524.010-031

NEEDLE-SHAPED DIAMOND DRILL BIT FOR REMOVING EPITHELIUM FROM HARD-TO-REACH PLACES, EVENING OUT THE SHAPE, ALSO IDEAL FOR COMBINED MANICURE.



172.514/524.012-040

DIAMOND DRILL BIT WITH THE SHAPE OF A FLAT END CONE FOR PUSHING, MATTIFYING, REMOVAL OF EPITHELIUM, FOR CORRECTING THE SHAPE.



257.514/524.012-025

THANKS TO ITS SMALL DIAMETER, IT ALLOWS YOU TO REMOVE CUTICLES EVEN DEEPER AND REMOVE NAIL MEMBRANES AT THE NAIL SHAFT MORE ACCURATELY. THE BIT IS PERFECT FOR SAFE LIFTING VERY DELICATE CUTICLES. IT IS INDISPENSABLE IN THE "PAINTING UNDER THE SKINS" TECHNIQUE.



243R.514/524.014-025

DIAMOND DRILL BIT IN THE SHAPE OF A LONG FLAME FOR: LIFTING THE CUTICLES, REMOVING THE EPITHELIA AND ALIGNING THE EDGES OF THE NAIL.



243XL.514/524.014-025

DIAMOND DRILL BIT IN THE SHAPE OF A LONG SHARP FLAME FOR: LIFTING THE CUTICLES, REMOVING THE EPITHELIA AND ALIGNING THE EDGES OF THE NAIL.



110.514/524.012-050

IT WILL SAFELY AND GENTLY MOVE AWAY THE CUTICLES AND GROWN EPITHELIA. IT GUARANTEES FAST AND PRECISE WORK.

New Expert drill bits



ZRN-266XUF HP.060

THE CHRISTMAS TREE-SHAPED CUTTER WITH CROSS-CUTS, WHICH WILL ENSURE THE SMOOTHNESS OF THE SURFACE, WILL ALSO BE AN IDEAL SOLUTION FOR BEGINNERS. THE UPPER PART OF THE CUTTER HAS VERY SMALL CUTS DESIGNED FOR SAFE WORK BY THE CUTICLES. DUE TO ITS DELICACY, IT WILL BE PERFECT FOR REMOVING GEL MANICURE.



ZRN-221RLM HP.040

IDEAL TOOL FOR CORRECTING DOWNWARD-GROWING NAILS AND ALSO FOR REMOVING NATURAL PLATE FOR GEL RENEWAL.

- THE BOLD HEAD PROTECTS THE PAD/OPIUM OF YOUR FINGER AND HYPONYHIUM FROM CUTTING AND UNPLEASANT BURNING.

THANKS TO SPECIAL LENGTH, THICKNESS AND BOLD HEAD ONE CAN MILL OUT ALL THE PLACES WITH DIFFICULT ACCESS.



THE CROSSCUTS ARE FACING BOTH SIDES, SO IT IS IDEAL FOR RIGHT- AND LEFT-HANDED PEOPLE.



ZRN-221RLF HP.060

THE CONE-SHAPED CARBIDE DRILL BIT WITH CUTS DEDICATED FOR BOTH RIGHT AND LEFT-HANDED STYLISTS. SAFE TOP THAT GUARANTEES SAFE WORK IN THE AREA OF THE NAIL FOLDS, EFFECTIVELY AND QUICKLY COPE WITH THE HYBRID AND GEL ON A NATURAL PLATE, HAS BEEN COVERED WITH A GOLD NANO-COATING.

Join IQNAILS world



MAILS@IQDENT.PL



IQDENT LTD. POLAND 93-124
ŁÓDŹ, KILINSKIEGO 228 STR.



WWW.IQNAILS.PL



@IQNAILSPOLAND



@IQNAILSPOLAND



@IQNAILSPOLAND



@IQNAILSPOLAND



@IQNAILSPOLAND



WWW.AKADEMIADIAMENTOWA.PL

